





ACAA1 Antibody

Product Code	CSB-PA001117LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09110
Immunogen	Recombinant Human 3-ketoacyl-CoA thiolase, peroxisomal protein (27-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-ketoacyl-CoA thiolase, peroxisomal (EC 2.3.1.16) (Acetyl-CoA acyltransferase) (Beta-ketothiolase) (Peroxisomal 3-oxoacyl-CoA thiolase), ACAA1, ACAA PTHIO
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ACAA1

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA001117LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA001117LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.