



LYZ Antibody, FITC conjugated

Product Code	CSB-PA02769C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61626
Immunogen	Recombinant Human Lysozyme C protein (19-148AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysozyme C (EC 3.2.1.17) (1,4-beta-N-acetylmuramidase C), LYZ, LZM
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LYZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.