





ISX Antibody

Product Code	CSB-PA011854LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2M1V0
Immunogen	Recombinant Human Intestine-specific homeobox protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intestine-specific homeobox (RAX-like homeobox), ISX, RAXLX
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ISX

Image



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA011854LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human glioma using CSB-PA011854LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.