



Polygalacturonase Antibody, HRP conjugated

Product Code	CSB-PA752116LB01CHAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7M1E7
Immunogen	Recombinant Chamaecyparis obtusa Polygalacturonase protein (51-514AA)
Raised In	Rabbit
Species Reactivity	Chamaecyparis obtusa
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polygalacturonase (PG) (EC 3.2.1.15) (Major pollen allergen Cha o 2) (Pectinase) (allergen Cha o 2)
Immunogen Species	Chamaecyparis obtusa
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.