





GGPS1 Antibody

Product Code	CSB-PA009393HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95749
Immunogen	Recombinant Human Geranylgeranyl pyrophosphate synthase protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Geranylgeranyl pyrophosphate synthase (GGPP synthase) (GGPPSase) (EC 2.5.1.-) ((2E,6E)-farnesyl diphosphate synthase) (Dimethylallyltranstransferase) (EC 2.5.1.1) (Farnesyl diphosphate synthase) (Farnesyltranstransferase) (EC 2.5.1.29) (Geranylgeranyl diphosphate synthase) (Geranyltranstransferase) (EC 2.5.1.10), GGPS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GGPS1

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA009393HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA009393HA01HU at dilution of 1:100