





RDH12 Antibody

Product Code	CSB-PA842727LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NR8
Immunogen	Recombinant Human Retinol dehydrogenase 12 protein (1-316AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinol dehydrogenase 12 (EC 1.1.1.-) (All-trans and 9-cis retinol dehydrogenase) (Short chain dehydrogenase/reductase family 7C member 2), RDH12, SDR7C2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RDH12

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA842727LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA842727LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.