



SNRPB2 Antibody, FITC conjugated

Product Code	CSB-PA362460LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08579
Immunogen	Recombinant Human U2 small nuclear ribonucleoprotein B" protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U2 small nuclear ribonucleoprotein B" (U2 snRNP B"), SNRPB2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SNRPB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.