





TRIM9 Antibody

Product Code	CSB-PA880153LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C026
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM9 protein (1-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM9 (EC 2.3.2.27) (RING finger protein 91) (RING-type E3 ubiquitin transferase TRIM9) (Tripartite motif-containing protein 9), TRIM9, KIAA0282 RNF91
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRIM9

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA880153LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA880153LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.