




NR1I2 Antibody

Product Code	CSB-PA016048ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75469
Immunogen	Recombinant Human Nuclear receptor subfamily 1 group I member 2 protein (175-434AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 1 group I member 2 (Orphan nuclear receptor PAR1) (Orphan nuclear receptor PXR) (Pregnane X receptor) (Steroid and xenobiotic receptor) (SXR), NR1I2, PXR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NR1I2

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA016048ESR1HU at dilution of 1:100