



entC3 Antibody, Biotin conjugated

Product Code	CSB-PA357923HD01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A0L5
Immunogen	Recombinant Staphylococcus aureus Enterotoxin type C-3 protein (28-266AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enterotoxin type C-3 (SEC3), entC3
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	entC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.