



ply Antibody, Biotin conjugated

Product Code	CSB-PA314690HD01FMW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C2J9
Immunogen	Recombinant Streptococcus pneumoniae serotype 4 Pneumolysin protein (2-471AA)
Raised In	Rabbit
Species Reactivity	Streptococcus pneumoniae serotype 4
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pneumolysin (Thiol-activated cytolysin), ply
Immunogen Species	Streptococcus pneumoniae serotype 4 (strain ATCC BAA-334 / TIGR4)
Research Area	Others
Target Names	ply
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.