



SCG5 Antibody

Product Code	CSB-PA020810DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05408
Immunogen	Recombinant Human Neuroendocrine protein 7B2 protein (27-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuroendocrine protein 7B2 (Pituitary polypeptide) (Secretogranin V) (Secretogranin-5) (Secretory granule endocrine protein I) [Cleaved into: N-terminal peptide; C-terminal peptide], SCG5, SGNE1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SCG5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.