





UCP1 Antibody

Product Code	CSB-PA025554ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25874
Immunogen	Recombinant Human Mitochondrial brown fat uncoupling protein 1 protein (1-307AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial brown fat uncoupling protein 1 (UCP 1) (Solute carrier family 25 member 7) (Thermogenin), UCP1, SLC25A7 UCP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UCP1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA025554ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA025554ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.