




CFHR3 Antibody

Product Code	CSB-PA005276ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02985
Immunogen	Recombinant Human Complement factor H-related protein 3 protein (80-200AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement factor H-related protein 3 (FHR-3) (DOWN16) (H factor-like protein 3), CFHR3, CFHL3 FHR3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CFHR3

Image



Western blot
 All lanes: Complement factor H-related protein 3 antibody at 6µg/ml + Rat heart tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 38, 31 kDa
 Observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.