




HSPB2 Antibody

Product Code	CSB-PA614956ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16082
Immunogen	Recombinant Human Heat shock protein beta-2 protein (1-182AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein beta-2 (HspB2) (DMPK-binding protein) (MKBP), HSPB2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPB2

Image




Western blot

All lanes: Heat shock protein beta-2 antibody at 3µg/ml + Mouse skeletal muscle tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 20 kDa
 Observed band size: 20 kDa



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA614956ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA614956ESR2HU at dilution of 1:100