



Minpp1 Antibody

Product Code	CSB-PA013831YA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O35217
Immunogen	Recombinant Rat Multiple inositol polyphosphate phosphatase 1 protein (31-481AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multiple inositol polyphosphate phosphatase 1 (EC 3.1.3.62) (2,3-bisphosphoglycerate 3-phosphatase) (2,3-BPG phosphatase) (EC 3.1.3.80) (Inositol (1,3,4,5)-tetrakisphosphate 3-phosphatase) (Ins(1,3,4,5)P(4) 3-phosphatase), Minpp1, Mipp1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Signal Transduction
Target Names	Minpp1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.