



LCP1 Antibody, Biotin conjugated

Product Code	CSB-PA012824YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13796
Immunogen	Recombinant Human Plastin-2 protein (2-627AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plastin-2 (L-plastin) (LC64P) (Lymphocyte cytosolic protein 1) (LCP-1), LCP1, PLS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LCP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.