





CBLIF Antibody

Product Code	CSB-PA009417YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27352
Immunogen	Recombinant Human Cobalamin binding intrinsic factor protein (19-417AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cobalamin binding intrinsic factor (Gastric intrinsic factor) (Intrinsic factor) (IF) (INF), CBLIF, GIF, IFMH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CBLIF


Image



Western blot
All lanes: CBLIF antibody at 2µg/ml + Mouse stomach tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 46, 43 kDa
Observed band size: 46 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA009417YA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA009417YA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.