





GRIA3 Antibody

Product Code	CSB-PA009900LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42263
Immunogen	Recombinant Human Glutamate receptor 3 protein (151-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor 3 (GluR-3) (AMPA-selective glutamate receptor 3) (GluR-C) (GluR-K3) (Glutamate receptor ionotropic, AMPA 3) (GluA3), GRIA3, GLUR3, GLURC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIA3

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA009900LA01HU at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA009900LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.