



GRIA3 Antibody, HRP conjugated

Product Code	CSB-PA009900LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42263
Immunogen	Recombinant Human Glutamate receptor 3 protein (151-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor 3 (GluR-3) (AMPA-selective glutamate receptor 3) (GluR-C) (GluR-K3) (Glutamate receptor ionotropic, AMPA 3) (GluA3), GRIA3, GLUR3, GLURC
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.