



PA Antibody, Biotin conjugated

Product Code	CSB-PA395880ED01ILR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4U6V9
Immunogen	Recombinant Influenza A virus Polymerase acidic protein (1-716AA)
Raised In	Rabbit
Species Reactivity	Influenza A virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polymerase acidic protein (EC 3.1.-.-) (RNA-directed RNA polymerase subunit P2), PA
Immunogen Species	Influenza A virus (strain A/USA:Huston/AA/1945 H1N1)
Research Area	Others
Target Names	PA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.