



# aroD Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA326162YB01SWW
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P24670
<b>Immunogen</b>	Recombinant Salmonella typhi 3-dehydroquinate dehydratase protein (1-252AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Salmonella typhi
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	3-dehydroquinate dehydratase (3-dehydroquinase) (EC 4.2.1.10) (Type I DHQase) (Type I dehydroquinase) (DHQ1), aroD
<b>Immunogen Species</b>	Salmonella typhi
<b>Research Area</b>	Others
<b>Target Names</b>	aroD
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.