



aroD Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA326162YD01SWW |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P24670 |
| Immunogen | Recombinant Salmonella typhi 3-dehydroquinase dehydratase protein (1-252AA) |
| Raised In | Rabbit |
| Species Reactivity | Salmonella typhi |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 3-dehydroquinase dehydratase (3-dehydroquinase) (EC 4.2.1.10) (Type I DHQase) (Type I dehydroquinase) (DHQ1), aroD |
| Immunogen Species | Salmonella typhi |
| Research Area | Others |
| Target Names | aroD |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |