





ALOXE3 Antibody

Product Code	CSB-PA863156ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYJ1
Immunogen	Recombinant Human Hydroperoxide isomerase ALOXE3 protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroperoxide isomerase ALOXE3 (EC 5.4.4.7) (Epidermis-type lipoxygenase 3) (Epidermal LOX-3) (e-LOX-3) (eLOX-3) (Hydroperoxy icosatetraenoate dehydratase) (EC 4.2.1.152), ALOXE3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ALOXE3

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA863156ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA863156ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.