



# dsbA Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA528735YA01FGP   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | O52376  |
| <b>Immunogen</b>           | Recombinant Pseudomonas syringae pv. tomato Thiol:disulfide interchange protein DsbA protein (23-214AA) |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Pseudomonas syringae pv. tomato   |
| <b>Tested Applications</b> | ELISA   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                        |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Thiol:disulfide interchange protein DsbA, dsbA  |
| <b>Immunogen Species</b>   | Pseudomonas syringae pv. tomato (strain ATCC BAA-871 / DC3000)  |
| <b>Research Area</b>       | Others  |
| <b>Target Names</b>        | dsbA  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                             |