



dsbA Antibody, Biotin conjugated

Product Code	CSB-PA528735YD01FGP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O52376
Immunogen	Recombinant Pseudomonas syringae pv. tomato Thiol:disulfide interchange protein DsbA protein (23-214AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas syringae pv. tomato
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thiol:disulfide interchange protein DsbA, dsbA
Immunogen Species	Pseudomonas syringae pv. tomato (strain ATCC BAA-871 / DC3000)
Research Area	Others
Target Names	dsbA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.