



entG Antibody, HRP conjugated

Product Code	CSB-PA362954HB01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A0L6
Immunogen	Recombinant Staphylococcus aureus Enterotoxin type G protein (26-258AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enterotoxin type G (SEG), entG, seg
Immunogen Species	Staphylococcus aureus (strain Mu50 / ATCC 700699)
Research Area	Others
Target Names	entG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.