



Prothoracicotropic hormone Antibody

Product Code	CSB-PA324609LA01BTT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17219
Immunogen	Recombinant Bombyx mori Prothoracicotropic hormone protein (116-224AA)
Raised In	Rabbit
Species Reactivity	Bombyx mori
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prothoracicotropic hormone (PTTH) [Cleaved into: P2K; P6K; Prothoracicotropic hormone]
Immunogen Species	BOMMO(Silk moth)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.