




KLRC2 Antibody

Product Code	CSB-PA012467LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26717
Immunogen	Recombinant Human NKG2-C type II integral membrane protein (94-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NKG2-C type II integral membrane protein (CD159 antigen-like family member C) (NK cell receptor C) (NKG2-C-activating NK receptor) (CD antigen CD159c), KLRC2, NKG2C
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KLRC2

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA012467LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.