



AQP4 Antibody, HRP conjugated

Product Code	CSB-PA11137B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55087
Immunogen	Recombinant Human Aquaporin-4 protein (253-323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-4 (AQP-4) (Mercurial-insensitive water channel) (MIWC) (WCH4), AQP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.