




CTSL Antibody

Product Code	CSB-PA006193LA01DSU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GKL8
Immunogen	Recombinant Chlorocebus aethiops Cathepsin L1 protein (114-333AA)
Raised In	Rabbit
Species Reactivity	Chlorocebus aethiops, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin L1 (EC 3.4.22.15) (Cathepsin L) (Major excreted protein) (MEP) [Cleaved into: Cathepsin L1 heavy chain; Cathepsin L1 light chain], CTSL, CTSL1
Immunogen Species	Chlorocebus aethiops
Research Area	Others
Target Names	CTSL

Image



Western blot

All lanes: CTSL antibody at 12μg/ml

Lane 1: A549 whole cell lysate

Lane 2: HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 38 kDa

Observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.