



X Antibody

Product Code	CSB-PA526046LA01HVP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O93195
Immunogen	Recombinant Hepatitis B virus genotype D Protein X protein (1-154AA)
Raised In	Rabbit
Species Reactivity	Hepatitis B virus genotype D
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein X (HBx) (Peptide X) (pX), X
Immunogen Species	Hepatitis B virus genotype D
Research Area	Others
Target Names	X