



Annexin A5 Antibody, Biotin conjugated

Product Code	CSB-PA001846HD01EQC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P70075
Immunogen	Recombinant Cynops pyrrhogaster Annexin A5 protein (1-323AA)
Raised In	Rabbit
Species Reactivity	Cynops pyrrhogaster
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A5 (Annexin V) (Annexin-5)
Immunogen Species	Cynops pyrrhogaster
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.