



entD Antibody, HRP conjugated

Product Code	CSB-PA324966LB01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20723
Immunogen	Recombinant Staphylococcus aureus Enterotoxin type D protein (26-258AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enterotoxin type D (SED), entD
Immunogen Species	Staphylococcus aureus
Research Area	Microbiology
Target Names	entD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.