



botB Antibody, HRP conjugated

Product Code	CSB-PA318083LB01CLQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10844
Immunogen	Recombinant Clostridium botulinum Botulinum neurotoxin type B protein (1-427AA)
Raised In	Rabbit
Species Reactivity	Clostridium botulinum
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Botulinum neurotoxin type B (BoNT/B) (EC 3.4.24.69) (Bontoxilysin-B) [Cleaved into: Botulinum neurotoxin B light chain; Botulinum neurotoxin B heavy chain], botB
Immunogen Species	Clostridium botulinum
Research Area	Others
Target Names	botB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.