



# botB Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA318083LD01CLQ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P10844
<b>Immunogen</b>	Recombinant Clostridium botulinum Botulinum neurotoxin type B protein (1-427AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Clostridium botulinum
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Botulinum neurotoxin type B (BoNT/B) (EC 3.4.24.69) (Bontoxilysin-B) [Cleaved into: Botulinum neurotoxin B light chain; Botulinum neurotoxin B heavy chain], botB
<b>Immunogen Species</b>	Clostridium botulinum
<b>Research Area</b>	Others
<b>Target Names</b>	botB
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.