




IL1A Antibody

Product Code	CSB-PA011613LA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48089
Immunogen	Recombinant Macaca mulatta Interleukin-1 alpha protein (113-271AA)
Raised In	Rabbit
Species Reactivity	Macaca mulatta, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 alpha (IL-1 alpha) (Hematopoietin-1), IL1A
Immunogen Species	Macaca mulatta
Research Area	Others
Target Names	IL1A

Image



Western blot
 All lanes: IL1A antibody at 6µg/ml
 Lane 1: Jurkat whole cell lysate
 Lane 2: Human placenta tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 31 kDa
 Observed band size: 31 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.