



TMEM27 Antibody, FITC conjugated

Product Code	CSB-PA881018EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBJ8
Immunogen	Recombinant Human Collectrin protein (15-141AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collectrin (Transmembrane protein 27), TMEM27
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMEM27
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.