



IL2RG Antibody, Biotin conjugated

Product Code	CSB-PA011651ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31785
Immunogen	Recombinant Human Cytokine receptor common subunit gamma protein (75-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytokine receptor common subunit gamma (Interleukin-2 receptor subunit gamma) (IL-2 receptor subunit gamma) (IL-2R subunit gamma) (IL-2RG) (gammaC) (p64) (CD antigen CD132), IL2RG
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL2RG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.