




ZAP70 (Ab-292) Antibody

Product Code	CSB-PA044449
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43403
Immunogen	Peptide sequence around aa.288-292 (N-S-D-G-Y) derived from Human ZAP70.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total ZAP70 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Clonality	Polyclonal
Alias	SRK, STD, TZK, ZAP-70
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZAP70

Image



Western blot analysis of extracts from 293, JK, p19 and K-562 cells using ZAP70(Ab-292) Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.