




# ZAP70 (Ab-292) Antibody

<b>Product Code</b>	CSB-PA044449
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P43403
<b>Immunogen</b>	Peptide sequence around aa.288-292 (N-S-D-G-Y) derived from Human ZAP70.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous level of total ZAP70 protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:1000
<b>Clonality</b>	Polyclonal
<b>Alias</b>	SRK, STD, TZK, ZAP-70
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ZAP70

## Image



Western blot analysis of extracts from 293, JK, p19 and K-562 cells using ZAP70(Ab-292) Antibody.