




MMP9 Antibody

Product Code	CSB-PA058909
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14780
Immunogen	Peptide sequence around aa.110-114(F-E-G-D-L) derived from Human MMP-9.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total MMP-9 protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Clonality	Polyclonal
Alias	GELB, CLG4B, MMP-9, MANDP2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MMP9

Image



Western blot analysis of extracts from HeLa cells using MMP-9 Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.