




TGFB1 Antibody

Product Code	CSB-PA286763
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01137
Immunogen	Peptide sequence around aa.378-382(Q-L-S-N-M) derived from Human TGF-b.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total TGF-b protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Clonality	Polyclonal
Alias	CED, LAP, DPD1, TGFB , TGFbeta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGFB1

Image



Western blot analysis of extracts from 293, HeLa and C2C12 cells using TGF-b Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.