




# MAPK11 (Ab-180/182) Antibody

<b>Product Code</b>	CSB-PA790629
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q15759
<b>Immunogen</b>	Peptide sequence around aa.179-183(M-T-G-Y-V) derived from Human P38 MAPK.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous level of total P38 MAPK protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:1000
<b>Clonality</b>	Polyclonal
<b>Alias</b>	RK, p38, CSBP, EXIP, Mxi2
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MAPK11

## Image



Western blot analysis of extracts from 293, 3T3 and PC12 cells using P38 MAPK(Ab-180/182) Antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.