




# CYCS Antibody

<b>Product Code</b>	CSB-PA877798
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P99999
<b>Immunogen</b>	Peptide sequence around aa.71-75(N-P-K-K-Y) derived from Human Cytochrome C.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous level of total Cytochrome C protein.
<b>Tested Applications</b>	ELISA,WB;WB:1:500-1:1000
<b>Clonality</b>	Polyclonal
<b>Alias</b>	CYC, HCS, THC4,
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CYCS

## Image



Western blot analysis of extracts from 3T3, C6 and A549 cells using Cytochrome C Antibody.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.