




MET (Ab-1234/1235) Antibody

Product Code	CSB-PA218313
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08581
Immunogen	Peptide sequence around aa.1232-1237(K-E-Y-Y-S-V) derived from Human Met.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total Met protein.
Tested Applications	ELISA,WB;WB:1:500-1:1000
Clonality	Polyclonal
Alias	HGF receptor, HGF-SF receptor, Met proto-oncogene tyrosine kinase, c-met, kinase Met
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MET

Image



Western blot analysis of extracts from HepG2 and HUVEC cells using Met (Ab-1234/1235) Antibody.