





Phospho-STAT5A (Tyr694) Antibody

Product Code	CSB-PA877372
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42229
Immunogen	Peptide sequence around phosphorylation site of tyrosine 694 (D-G-Y(p)-V-K) derived from Human STAT5A.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of STAT5A only when phosphorylated at tyrosine 694.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:100
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	MGF; MPF; Mammary gland factor; STA5A,; STAT5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	STAT5A

Image



Western blot analysis of extracts from K562 cells untreated or treated with G-CSF using STAT5A(Phospho-Tyr694) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using STAT5A(Phospho-Tyr694) Antibody(left) or the same antibody preincubated with blocking peptide(right).

Product Modify Phospho-Tyr694



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.