





Phospho-EZR (Tyr353) Antibody

Product Code	CSB-PA988037
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15311
Immunogen	Peptide sequence around phosphorylation site of tyrosine 353 (Q-D-Y(p)-E-E) derived from Human EZRIN.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of Ezrin only when phosphorylated at tyrosine 353.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:100
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography
Clonality	Polyclonal
Alias	Cytovillin; EZRI; VIL2; Villin 2; p81
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EZR

Image



Western blot analysis of extracts from 293 cells untreated or treated with EGF using Ezrin(Phospho-Tyr353) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Ezrin(Phospho-Tyr353) Antibody(left) or the same antibody preincubated with blocking peptide(right).

Product Modify Phospho-Tyr353