





Phospho-PTK2 (Tyr925) Antibody

Product Code	CSB-PA262237
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05397
Immunogen	Peptide sequence around phosphorylation site of tyrosine 925 (K-V-Y(p)-E-N) derived from Human FAK.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of FAK only when phosphorylated at tyrosine 925.
Tested Applications	ELISA,WB,IF;WB:1:500-1:1000,IF:1:100-1:200
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	FADK 1; FAK1; PTK2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTK2

Image



Western blot analysis of extracts from 293 cells using FAK(Phospho-Tyr925) Antibody(Lane 2) and the same antibody preincubated with blocking peptide(Lane1).



Immunofluorescence staining of methanol-fixed Hela cells using FAK(Phospho-Tyr925) Antibody.

Product Modify

Phospho-Tyr925

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.