





Phospho-RB1 (Ser807) Antibody

Product Code	CSB-PA193785
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06400
Immunogen	Peptide sequence around phosphorylation site of serine 807 (Y-I-S(p)-P-L) derived from Human Rb.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of Rb only when phosphorylated at serine 807.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:1000,IHC:1:50-1:100
Form	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography usi
Clonality	Polyclonal
Alias	P105-RB; PP105; PP110; RB-1; RB1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RB1

Image



Western blot analysis of extracts from 3T3 cells untreated or treated with UV using Rb(Phospho-Ser807) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Rb(Phospho-Ser807) Antibody(left) or the same antibody preincubated with blocking peptide(right).



Product Modify

Phospho-Ser807

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.